



CooperVision®  
Live Brightly.®

# Product Specifications Guide



## MyDay® daily disposable family

### MyDay® daily disposable



#### Wearing Schedule

Daily disposable

#### Carton Size Options

90 pack

#### Material / H<sub>2</sub>O Content

stenfilcon A / 54%

**Base Curve** 8.4mm

**Diameter** 14.2mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** Yes

#### Sphere Power

+8.00D to -12.00D  
(0.50D steps after  
+5.00D and  
-6.00D)

#### Design

Aspheric

High-performing 1 Day contact lenses for demanding lifestyles.

### MyDay® toric daily disposable



#### Wearing Schedule

Daily disposable

#### Carton Size Options

90 pack

#### Material / H<sub>2</sub>O Content

stenfilcon A / 54%

**Base Curve** 8.6mm

**Diameter** 14.5mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** Yes

#### Sphere Power

Plano to -6.00D  
(0.25 steps)

#### Cylinder Power

-0.75  
-1.25  
-1.75

#### Axis

10°, 20°, 30°, 40°,  
50°, 60°, 70°, 80°,  
90°, 100°, 110°,  
120°, 130°, 140°,  
150°, 160°, 170°, 180°

#### Design

Optimised  
Toric Lens  
Geometry™

Plano to -6.00D  
(0.25D steps)

-2.25

10°, 20°, 90°,  
160°, 170°, 180°

-6.50 to -10.00D  
(0.50D steps)

+0.50D to +6.00D  
(0.50 steps)

-6.50D to -10.00D  
(0.50D steps)+  
0.25D to +6.00D  
(0.25D steps)

-0.75  
-1.25  
-1.75

10°, 20°, 70°, 80°,  
90°, 100°, 110°,  
160°, 170°, 180°

High-performing 1-day contact lenses for wearers with astigmatism.

## clariti® 1 day family

### clariti® 1 day



#### Wearing Schedule

Daily disposable

#### Carton Size Options

90 pack

#### Material / H<sub>2</sub>O Content

somofilcon A / 56%

**Base Curve** 8.6mm

**Diameter** 14.1mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** No

#### Sphere Power

-0.50D to -10.00D  
+0.50D to +8.00D  
(0.50D steps after  
+/-6.00D)

#### Design

Aspheric

High breathability<sup>1</sup>, comfort and great value<sup>†</sup> in convenient daily disposable contact lenses.

### clariti® 1 day toric



#### Wearing Schedule

Daily disposable

#### Carton Size Options

30 pack

#### Material / H<sub>2</sub>O Content

somofilcon A / 56%

**Base Curve** 8.6mm

**Diameter** 14.3mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** No

#### Sphere Power

Plano to -6.00D  
(0.25D steps)

#### Cylinder Power

-0.75  
-1.25  
-1.75

#### Axis

10° to 180°  
(10° steps)

#### Design

Aspheric

High breathability<sup>1</sup>, comfort and great value<sup>†</sup> in convenient daily disposable contact lenses.

Plano to -6.00D  
(0.25D steps)

-2.25

10°, 20°, 70°,  
80°, 90°, 100°,  
110°, 160°, 170°,  
180°

+0.25D to +4.00D  
(0.25D steps)

-0.75  
-1.25  
-1.75

10°, 20°, 70°,  
80°, 90°, 100°,  
110°, 160°, 170°,  
180°

-6.50 to -9.00  
(0.50D steps)

-0.75  
-1.25  
-1.75

10°, 20°, 60°,  
70°, 80°, 90°,  
100°, 110°, 120°,  
160°, 170°, 180°

-6.50 to -9.00  
(0.50D steps)

-2.25

10°, 20°, 90°,  
160°, 170°, 180°

## clariti® 1 day multifocal



High breathability<sup>1</sup>, comfort and great value<sup>†</sup> in convenient daily disposable contact lenses.

**Wearing Schedule**  
Daily disposable

**Sphere Power**  
+5.00D to -6.00D  
(0.25D steps)

**Carton Size Options**  
30 pack

**Add Power**  
Low:  
up to +2.25D

High:  
+2.25D to  
+3.00D

**Material / H<sub>2</sub>O Content**  
somofilcon A / 56%

**Design**  
Simultaneous vision

**Base Curve** 8.6mm  
**Diameter** 14.1mm  
**FDA Group** SiHy  
**UV Blocking\*** Yes  
**Visibility Tint** No

## Biofinity® family

### Biofinity Energys®



The world's only contact lenses with Digital Zone Optics™. Designed for today's digital life.

**Wearing Schedule**  
Daily wear or extended wear up to 29 nights/30 days

**Sphere Power**  
+8.00D to -12.00D  
(0.50D steps after +/-6.00D)

**Carton Size Options**  
6 pack

**Design**  
Aspheric  
Digital Zone Optics® lens design

**Material / H<sub>2</sub>O Content**  
comfilcon A / 48%

**Base Curve** 8.6mm  
**Diameter** 14.0mm  
**FDA Group** SiHy  
**UV Blocking\*** No  
**Visibility Tint** Yes

### Biofinity®



Premium monthly replacement contact lenses for long-lasting comfort and excellent vision<sup>2-4</sup>

**Wearing Schedule**  
Daily wear or extended wear up to 29 nights/30 days

**Sphere Power**  
+8.00D to -12.00D  
(0.50D steps after +/-6.00D)

+8.50D to +15.00D  
-12.50D to -20.00D  
(0.50D steps)

**Carton Size Options**  
6 pack

**Design**  
Aspheric

**Material / H<sub>2</sub>O Content**  
comfilcon A / 48%

**Base Curve** 8.6mm  
**Diameter** 14.0mm  
**FDA Group** SiHy  
**UV Blocking\*** No  
**Visibility Tint** Yes

## Biofinity® toric



The perfect solution for astigmatism. With Biofinity® toric monthly contact lenses you can enjoy excellent vision and a more comfortable lens-wearing experience.<sup>2,5-6</sup>

### Wearing Schedule

Daily wear or extended wear up to 29 nights/30 days

### Carton Size Options

6 pack

### Material / H<sub>2</sub>O Content

comfilcon A / 48%

**Base Curve** 8.7mm

**Diameter** 14.5mm

**FDA Group** SiHy

**UV Blocking\*** No

**Visibility Tint** Yes

### Design

Optimised  
Toric Lens  
Geometry™

### Sphere Power

+8.00D to -10.00D  
(0.50D steps after  
+/-6.00D)

### Cylinder Power

-0.75  
-1.25  
-1.75  
-2.25

### Axis

10° to 180°  
(10° steps)

+20.00D to -20.00D  
(0.50D steps after  
+/-6.00D)

-2.75 to -5.75  
(0.50D steps)

5° to 180°  
(5° steps)

+8.50 to +20.00  
(0.50D steps)

-0.75  
-1.25  
-1.75  
-2.25

5° to 180°  
(5° steps)

-10.50 to -20.00  
(0.50D steps)

## Biofinity® multifocal



Premium monthly replacement contact lenses for presbyopia, to provide you incredible comfort<sup>1</sup> and excellent vision at all distances<sup>7</sup>

### Wearing Schedule

Daily wear or extended wear up to 29 nights/30 days

### Carton Size Options

6 pack

### Material / H<sub>2</sub>O Content

comfilcon A / 48%

**Base Curve** 8.6mm

**Diameter** 14.0mm

**FDA Group** SiHy

**UV Blocking\*** No

**Visibility Tint** Yes

### Sphere Power

+6.00D to -10.00D  
(0.50D steps after  
-6.00D)

### Add Power

+1.00  
+1.50  
+2.00  
+2.50

### Design

D lens  
N lens

## Biofinity® toric multifocal



Designed to keep up with your evolving vision needs.

### Wearing Schedule

Daily wear for monthly replacement

### Carton Size Options

6 pack

### Material / H<sub>2</sub>O Content

comfilcon A / 48%

**Base Curve** 8.7mm

**Diameter** 14.5mm

**FDA Group** SiHy

**UV Blocking\*** No

**Visibility Tint** Yes

### Sphere Power

+10.00 to -10.00  
(0.50D steps after  
+/-6.00)

### Cylinder Power

-0.75 to -5.75  
(0.50D steps)

### Axis


5° to 180°  
(5° steps)

### Add Power

+1.00 to  
+2.50  
(0.50D steps)

## Serenity™ family

### Serenity™



High performance. Outstanding value.

#### Wearing Schedule

Daily

#### Sphere Power

+8.00D to -12.00D  
(0.50D steps after  
+/-6.00D)

No plano

#### Carton Size Options

6 pack

#### Design

Aspheric front surface

#### Material / H<sub>2</sub>O Content

fanfilcon A / 55%

**Base Curve** 8.4mm


**Diameter** 14.2mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** Yes

### Serenity™ toric



High performance. Outstanding value.

#### Wearing Schedule

Daily

#### Sphere Power

-10.00D to +8.00D  
(0.50D steps after  
+/- 6.00)

#### Carton Size Options

6 pack

#### Cylinder Power

-0.75 to -2.25

#### Material / H<sub>2</sub>O Content

fanfilcon A / 55%

#### Axis

10° to 180°  
(10° steps)

#### Design

Optimised Toric  
Lens Geometry™

**Base Curve** 8.5mm

**Diameter** 14.5mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** Yes

## Brilliant Futures™ with MiSight® Myopia Management

### Brilliant Futures™ with MiSight® Myopia Management



Child-friendly daily disposable contact lenses proven to slow the progression of short-sightedness in children.<sup>8</sup>

#### Wearing Schedule

Daily

#### Sphere Power

-6.00D to -0.50D  
(0.25D steps)

#### Carton Size Options

90 pack

#### Material / H<sub>2</sub>O Content

omafilcon A

**Base Curve** 8.7mm

**Diameter** 14.2mm

**FDA Group** 2

**UV Blocking\*** No

**Visibility Tint** Yes

## Proclear® family

### Proclear® 1 day



Clear vision<sup>†</sup>, convenience and comfort, for contact lens wearers who experience eye dryness.

<b>Wearing Schedule</b>	<b>CARTON SIZE OPTIONS</b>	<b>Material / H<sub>2</sub>O Content</b>	<b>Base Curve</b>	8.7mm
Daily	30 pack 90 pack	omafilcon A / 60%	<b>Diameter</b>	14.2mm
<b>Sphere Power</b>	<b>Design</b>		<b>FDA Group</b>	2
+8.00D to -12.00D (0.50D steps after -6.00D and +5.00D)	Aspheric		<b>UV Blocking*</b>	No
			<b>Visibility Tint</b>	Yes

### Proclear® 1 day multifocal



Convenience and comfort for multifocal contact lens wearers who experience eye dryness.

<b>Wearing Schedule</b>	<b>Carton Size Options</b>	<b>Material / H<sub>2</sub>O Content</b>	<b>Base Curve</b>	8.7mm
Daily disposable	30 pack 90 pack	omafilcon A / 60%	<b>Diameter</b>	14.2mm
<b>Sphere Power</b>	<b>Add Power</b>	<b>Design</b>	<b>FDA Group</b>	2
+6.00D to -10.00D (0.50D steps after -6.00D)	Single power profile; accommodates up to +2.50 ADD	Centre near aspheric	<b>UV Blocking*</b>	No
			<b>Visibility Tint</b>	Yes

### Proclear®



Clear vision<sup>9</sup> and comfort for contact lens wearers who experience eye dryness.

<b>Wearing Schedule</b>	<b>Carton Size Options</b>	<b>Material / H<sub>2</sub>O Content</b>	<b>Base Curve</b>	8.2, 8.6mm
Daily	6 pack	omafilcon B / 62%	<b>Diameter</b>	14.2mm
<b>Sphere Power</b>	<b>Design</b>		<b>FDA Group</b>	2
+20.00D to -20.00D (0.50D steps after +/-6.00D)	Sphere		<b>UV Blocking*</b>	No
Plano to -10.00D (0.50D steps after -6.00D)			<b>Visibility Tint</b>	Yes

## Proclear® toric



Helps provide all-day comfort for contact lens wearers with astigmatism who experience eye dryness.

<b>Wearing Schedule</b> Daily	<b>Carton Size Options</b> 6 pack	<b>Material / H<sub>2</sub>O Content</b> omafilcon B / 62%	<b>Base Curve</b> 8.4, 8.8mm <b>Diameter</b> 14.4mm <b>FDA Group</b> 2 <b>UV Blocking*</b> No <b>Visibility Tint</b> Yes
<b>Sphere Power</b> +6.00D to -8.00D (0.50D steps after -6.50D)	<b>Cylinder Power</b> -0.75 -1.25 -1.75 -2.25	<b>Axis</b> 10° to 180° (10° steps)	
+10.00D to -10.00D (0.50D steps after +/-6.50D)	-0.75 to -5.75 (0.50D steps)	5° to 180° (5° steps)	

## Proclear® multifocal



Helps provide all-day comfort for multifocal contact lens wearers who experience eye dryness.

<b>Wearing Schedule</b> Daily	<b>Carton Size Options</b> 6 pack	<b>Material / H<sub>2</sub>O Content</b> omafilcon B / 62%	<b>Base Curve</b> 8.4, 8.7mm <b>Diameter</b> 14.4mm <b>FDA Group</b> 2 <b>UV Blocking*</b> No <b>Visibility Tint</b> Yes
<b>Sphere Power</b> +6.00D to -8.00D (0.50D steps after -6.50D)	<b>Add Power</b> +1.00 +1.50 +2.00 +2.50	<b>Design</b> D lens N lens	
+20.00D to -20.00D (0.50D steps after +/-6.50D)	+1.00 to +4.00 (0.50D steps)		

## Proclear® multifocal toric



Helps provide all-day comfort for multifocal contact lens wearers with astigmatism who experience eye dryness.

<b>Wearing Schedule</b> Daily	<b>Carton Size Options</b> 6 pack	<b>Material / H<sub>2</sub>O Content</b> omafilcon B / 62%	<b>Base Curve</b> 8.4, 8.8mm <b>Diameter</b> 14.4mm <b>FDA Group</b> 2 <b>UV Blocking*</b> No <b>Visibility Tint</b> Yes
<b>Sphere Power</b> +20.00D to -20.00D (0.50D steps after +/-6.50D)	<b>Cylinder Power</b> -0.75 to -5.75 (0.50D steps)	<b>Axis</b> 5° to 180° (5° steps)	
<b>Design</b> D lens N lens			



## Biomedics<sup>®</sup> family

### Biomedics<sup>®</sup> 55 Evolution



Designed for comfort and performance at an affordable price.

**Wearing Schedule**  
Daily

**Carton Size Options**  
6 pack

**Material / H<sub>2</sub>O Content**  
ocufilcon D / 55%


**Base Curve** 8.6, 8.8, 8.9mm  
**Diameter** 14.2, 14.4mm  
**FDA Group** 4  
**UV Blocking\*** Yes  
**Visibility Tint** Yes

**Sphere Power**  
-0.25D to -10.00D  
(0.50D steps after -6.00D)

**Design**  
Aspheric

+0.25D to +6.00D  
(0.50D steps after +5.00D)

### Biomedics<sup>®</sup> toric



The affordable contact lens for wearers with astigmatism.

**Wearing Schedule**  
Daily

**Carton Size Options**  
6 pack

**Material / H<sub>2</sub>O Content**  
ocufilcon D / 55%

**Base Curve** 8.7mm  
**Diameter** 14.5mm  
**FDA Group** 4  
**UV Blocking\*** Yes  
**Visibility Tint** Yes

**Sphere Power**  
+6.00D to -9.00D  
includes Plano  
(0.50D steps after -6.00D and +5.00D)

**Cylinder Power**  
-0.75  
-1.25  
-1.75  
-2.25

**Axis**  
10° to 180°  
(10° steps)

<sup>†</sup> During daily wear..

<sup>‡</sup> Due to aspheric optics.

<sup>1</sup> Brennan NA. Beyond Flux: Total Corneal Oxygen Consumption as an Index of Corneal Oxygen During Contact Lens Wear. Optom Vis Sci 2005. 2005;82(6):467-472.

<sup>2</sup> CooperVision data on file. 2006. Aquaform<sup>®</sup> Technology uses a matrix of long silicone chains with hydrogen bonds to lock water molecules in the lens for incredible, long-lasting comfort and clarity.

<sup>3</sup> CooperVision data on file. 2014. Clinical evaluation of Biofinity<sup>®</sup> spherical contact lens compared to Air Optix<sup>®</sup> Aqua spherical lenses over one month of wear.

<sup>4</sup> CooperVision data on file. 2014. Clinical evaluation of Biofinity<sup>®</sup> spherical contact lens compared to Ultra spherical contact lens over one month of wear.

<sup>5</sup> CooperVision data on file. 2014. Clinical evaluation of Biofinity<sup>®</sup> toric compared to AIR OPTIX<sup>®</sup> for ASTIGMATISM contact lens over one month of wear.

<sup>6</sup> CooperVision data on file, 2016. Clinical evaluation of Biofinity<sup>®</sup> toric compared to ULTRA<sup>®</sup> for astigmatism over one month of wear.

<sup>7</sup> CooperVision data on file 2011. Clinical Evaluation of Biofinity multifocal<sup>®</sup> compared to ACUVUE<sup>®</sup> Oasys for Presbyopia, AIR OPTIX<sup>®</sup> for Presbyopia and PureVision<sup>®</sup> Multi-Focal.

<sup>8</sup> Chamberlain P, et al. A 3-year randomized clinical trial of MiSight<sup>®</sup> lenses for myopia control. Optom Vis Sci. 2019;96(8):556-567.

<sup>9</sup> CCVI data on file 2017 & 2014; non-dispensing short-term fitting study. N=40 & n=20.

\*UV-blocking contact lenses help provide protection against transmission of harmful UV radiation to eye but are not substitutes for protective UV-absorbing eyewear, such as UV absorbing goggles or sunglasses, as they do not completely cover the eye or surrounding area. Continue to use UV-absorbing eyewear as directed by your eye care professional.